**RS232 Interface**  
INSTRUCTION MANUAL

**Interfaz RS232**  
MANUAL DE INSTRUCCIONES

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## INTRODUCTION


This Interface Kit is for use with OHAUS Scout STX, SPX, SKX, SJX series and PJX series products.

## INTERFACE INSTALLATION

Install the Interface module on the mini DIN port (round connector) located at the rear side of the balance as shown. During this process, please make sure the balance is powered off.

### Notes:

- Please make sure the small triangle mark (▲) on the round connector is at the bottom and then plug in.
- The appearance of your model may be different.

## SETUP

- A) The OHAUS Interface is preset to communicate using the following settings: **9600 baud, 8 bit, no parity, no handshake**. Determine the RS232 parameters required for the printer or computer that is to be connected. See computer or printer documentation if assistance is required. If the parameters do not match, it will be necessary to change either the balance settings, or the computer / printer settings.
- B) Upon installation of the interface, the balance will recognize the RS232 Interface and add relevant items to the menu. Configure the balance to the desired RS232 and printing parameters; refer to the balance Instruction Manual for assistance in using the menus.

**For SPX, SKX and SJX Balances:**

Upon installation, the balance will show [r.5.2.3.z] during the power sequence.

<p><b>PRINT</b></p> <p>Reset: <b>no</b>, yes</p> <p>Stable: <b>off</b>, on</p> <p>A.Print: <b>off</b></p> <p>On.Stable interval accept continuous</p> <p>Content: Result (-&gt; <b>off</b>, on) Gross (-&gt; <b>off</b>, on) Net (-&gt; <b>off</b>, on) Tare (-&gt; <b>off</b>, on) Header (-&gt; <b>off</b>, on) Footer (-&gt; <b>off</b>, on) Mode (-&gt; <b>off</b>, on) Unit (-&gt; off, <b>on</b>) Info (-&gt; <b>off</b>, on)</p> <p>Layout: Format (-&gt;Single,<b>Multi</b>)</p> <p>Feed (-&gt;<b>Line</b>, 4 Lines, Form)</p> <p>End Print: Exit menu</p>	<p><b>RS232</b></p> <p>Reset: <b>no</b>, yes</p> <p>Baud Rate: 1200...<b>9600</b>...115200</p> <p>Parity: 7 even, 7 odd, 7 none, <b>8 none</b></p> <p>Stop bit: <b>1</b>, 2</p> <p>Handshake: <b>none</b>, XOn-XOff</p> <p>Alternate command: Print (A...<b>P</b>...Z) Tare (A...<b>T</b>...Z) Zero (A...<b>Z</b>)</p> <p>End RS232: Exit menu</p>
--	--




**Note:** Items with bold font are default settings.

**SETUP DEFINITIONS**

Menu Settings	Explanation
PRINT / A. Print – Off	Only sends data when PRINT is pressed.
PRINT / A.Print / On.Stable – Load, Load.Zero	Balance only transmits stable data. <i>Load:</i> Print stable value excluding zero <i>Load.Zero:</i> Print stable value including zero
PRINT / A.Print / interval - (xx) sec	Balance sends data every (xx) seconds.
PRINT / Auto Print – accept	Balance only sends stable accept data in Check mode.
PRINT / A.Print – Continuous	Balance repeatedly sends data as fast as possible.
PRINT / Content	On: print relevant content after weighing data Gross – on: G; Net – on: N; Tare – on: T; Info – on: balance will print application mode setting parameters
PRINT / Layout / Format - (Single, Multi)	<i>Single:</i> print all data in one line <i>Multi:</i> print all data in multiple lines
PRINT / Layout / Feed - (Line, 4 Lines, Form)	<i>Line:</i> feed one line after printing <i>4 Lines:</i> feed four lines after printing <i>Form:</i> feed one page after printing (move to the top of next page after printing)

**For STX Balances:**

Upon installation, the STX balance will show RS232 icon () in the upper right corner.

Menu added	Function
Communication  --RS232  --Baud Rate	Set the baud rate (bits per second). 
--Transmission	Set the data bits, stop bit, and parity. 
--Handshake	Set the flow control method. 
--Print Settings	For more info, please refer to the STX instruction manual. 

**For PJX Balances:**


Menu Navigation	
<p><b>PRINT</b>                      Header                      Date and Time                      Balance ID                      Balance Name                      User Name                      Project Name                      Application Name                      Result                      Gross Weight                      Net Weight                      Tare Weight                      Line Feed</p>	<p><b>RS232</b>                      Baud Rate                      Transmission                      Handshake</p>

Menu Settings	
<p>RS232 Interface Setup</p> <p>--Baud Rate</p> <p>--Transmission</p> <p>--Handshake</p>	<p>Enter this sub-menu to customize RS232 standard settings. Data may be output to either a printer or a PC.</p> <p>Set the baud rate (bits per second).</p> <p>1200 = 1200 bps                      2400 = 2400 bps                      4800 = 4800 bps  <b>9600 = 9600 bps</b>                      19200 = 19200 bps                      38400 = 38400 bps</p> <p>Set the data bits, stop bit, and parity.</p> <p><b>8-N-1 = 8 data bits, no parity, stop bit 1</b>                      8-N-2 = 8 data bits, no parity, stop bit 2                      7-E-1 = 7 data bits, even parity, stop bit 1                      7-E-2 = 7 data bits, even parity, stop bit 2                      7-N-1 = 7 data bits, no parity, stop bit 1                      7-N-2 = 7 data bits, no parity, stop bit 2                      7-O-1 = 7 data bits, odd parity, stop bit 1                      7-O-2 = 7 data bits, odd parity, stop bit 2</p> <p>Set the flow control method.</p> <p><b>NONE = no handshaking</b>                      XON-XOFF = XON/XOFF handshaking                      HARDWARE = hardware handshaking</p>
<p>Print Settings</p>	<p>For more info, please refer to the PJX instruction manual.</p>

## RS232 (DB9) Pin Connections

The RS232 Interface Kit cable terminates with a 9 pin Sub-D female connector. Active pins: Pin 2 = TXD, Pin 3 = RXD, Pin 5 = Ground.

The Interface connector can connect directly to a computer, an OHAUS printer and other printers.

Diagram	Type	Description
 <p>Note: RTS, CTS and RT_CTL signal pins are only available for PJX balances.</p>	Interface type	Voltage interface conforming to EIA RS-232C/DIN 66020 (CCITT V24/V.28)
	Max. cable length	15 m
	Signal level	Output: +5 V ... +15 V (RL = 3 – 7kΩ) -5 V ... -15 V (RL = 3 - 7 kΩ) Input: +3 V ... +25 V -3 V ... -25 V
	Connector	Sub-D, 9-pole, female
	Operating mode	Full duplex
	Transmission mode	Bit-serial, asynchronous
	Transmission code	ASCII
	Baud rates	Scout STX, SPX, SKX, SJX balances: 1200, 2400, 4800, 9600, 19200, 38400, 57600, 115200 (firmware selectable)
	Bits/parity	7-bit/even, 7-bit/odd, 7-bit/none, 8-bit/none (firmware selectable)
	Stop bits	Stop bit 1, 2
	Handshake	None, XON/XOFF, RTS/CTS (selectable)
	End-of-line	Not selectable

## RS232 OUTPUT

With the interface installed the balance will operate in several ways according to the settings.

Print format can be switched by xFMT user command (see command table in RS232 INPUT section).

### New Scout Print Format (Default Format):

#### Output String (Non Check Weighing Applications):

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[T/N/G/PT]	2 characters (right justified)
[Term]	2 characters

**Note:** All of the fields have fixed length.

#### Output example:

```

*****192.21_****g_*_*
*****0.01_****g_?_*

*****95.0_****g_*_*N
*****169.6_****g_*_*G
*****95.0_****g_*_*N
*****74.6_****g_*_*T

```

A.Print: off; Stable: off  
Stable reading  
unstable reading

A.Print: off; Stable: on  
Content / Result ->  
Content / Gross ->  
Content / Net ->  
Content / Tare ->

#### Output String (Check Weighing Application):

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[T/N/G/PT]	2 characters (right justified)
[space]	1 character
[application status]	6 characters (right justified)
[Term]	2 characters

Example:

```
*****192.21_****g_***_Accept
*****0.01_****g_?_*_Under
```

A.Print: off; Stable: off  
 Stable reading,  
 Unstable reading

**Scout Pro Print Format 1 (for the models 303/123/202/402/602/2001/6001/401FZH/601FZH/6000FZH):**

Output String:

[weight]	12 characters (right justified)
[space]	1 character
[unit]	5 characters (left justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[Legend]	1~10 characters

Output example:

```
*****0.00_g****_*
*****12.73_g****_?
*****0.85_oz***_WET*WT
```

**Scout Pro Print Format 2 (for the models 401/601/6000):**

Output String:

[weight]	11 or 12 characters (right justified)
[space]	1 character
[unit]	1~5 characters
[space]	1 character
[stability indicator]	1 character; "?" when unstable, space when stable
[space]	1 character
[Legend]	1~10 characters

**Note:** The unit field length varies with different units. The weight field could be 11 or 12, depending on if the weight string has a dot or not.

Output example:

```
*****100_g_*
*****273_g_?_*
*****8.5_oz***_WET*WT
```

**Print Format 3 (for Certain POS Systems):**

Output String:

[weight]	11 characters (right justified)
[space]	1 character
[unit]	5 characters (right justified)
[stability indicator]	1 character; "?" when unstable, space when stable
[Term]	2 characters

**Note:** All of the fields have fixed length.

Output example:

```

*****0.00_****g*
*****12.73_****g?
    
```

**PJX Print Format:**

The Result Data, and G/N/T data, is output in the following format.

Field:	Label	Space	Weight	Space	Unit	Space	Stability	Space	G/N	Space	Term. Characters
Length:		1	11	1	5	1	≤ 1	≤ 1	≤ 3	0	≤ 8

## RS232 INPUT

The following interface commands will be acknowledged by the Balance. They are case sensitive. The balance will return “ES” for invalid commands.

For SPX, SKX, SJX and STX Balances:

Command	Function
IP	Immediate Print of displayed weight (stable or unstable).
P	Print displayed weight (stable or unstable).
CP	Continuous Print.
SP	Print on Stability.
SLP	Auto Print stable non-zero displayed weight.
SLZP	Auto Print stable non-zero weight and stable zero reading.
xP	Interval Print x = Print Interval (1-3600 sec) 0P ends interval Print
0P	Turn off Auto Print
H	<b>H x “text”</b> Enter Header line , where x = line number 1 to 5, “text” = header text up to 24 alphanumeric characters
F	<b>F x “text”</b> Enter Footer line , where x = line number 1 to 2, “text” = footer text up to 24 alphanumeric characters
Z	Same as pressing Zero Key
T	Same as pressing Tare Key.
xT	Establish a preset Tare value in displayed unit. x = preset tare value. Sending 0T clears tare (if allowed).
PT	Prints Tare weight stored in memory.
PM	Print current application mode (weighing mode).
xM	Set current application mode to x. x depends on application 1M: WEIGHT, 2M: COUNT, 3M: PERCENT, 4M: CHECK, 5M: DYNAMIC, 6M: TOTAL, 7M: DENSITY, 8M: HOLD, 9M: MOLE
M	Scroll to the next enabled mode.
PU	Print Current weighing unit: g, kg, lb, oz, etc....
xU	Set balance to unit x: g, kg etc. 1U: g, 2U: kg, 3U: ct, 4U: N, 5U: oz, 6U: ozt, 7U: dwt, 8U: lb, 9U: lb:oz, 10U: gm, 11U: thk, 12U: tsg, 13U: ttw , 14U: tola, 15U: c
U	Scroll to the next enabled unit.
ON	Brings out of Standby
OFF	Goes to Standby.
C	Begin Span Calibration
AC	Abort Calibration.
PSN	Print Serial Number.
PV	Print Version: print name, software revision and LFT ON (if LFT is set ON).
x#	Set Counting APW (x) in grams. (must have APW stored)
P#	Print Counting application APW.
x%	Set Percent application reference weight (x) in grams. (must have reference weight stored)
P%	Print Percent application reference weight.
xCO	Set Checkweighing Over Limit in grams x.
xCU	Set Checkweighing Under Limit in grams x.
PCO	Print Checkweighing Over Limit.
PCU	Print Checkweighing Under Limit.
xMM	Set Molar Mass in g/mol .
PMM	Print Molar Mass
xS	0 = print unstable data, 1 = print stable only
xFMT	0 = New Scout print format (default) ; 1 = Scout Pro print format 1; 2 = Scout Pro print format 2; 3 = for certain POS system.
xRL	0 = disable response; 1 = enable response.

## For PJX balances:

Command	Function
IP	Immediate Print of displayed weight (stable or unstable).
P	Print displayed weight (stable or unstable).
CP	Continuous Print.
SP	Print on Stability.
SLP	Set Auto Print to On Stability, allow non-zero displayed weight be printed.
SLZP	Set Auto Print to On Stability, allow both stable non-zero weight and stable zero reading printed.
xP	Set Auto Print to Interval Print, x = print interval (1-3600 sec), 0P disable the interval Print
0P	0P disable interval print, continuous print or print on stability
H	Enter Print Header Lines, the format is: H x "header string". Where x = line number 1 to 3, "header string" can be up to 24 alphanumeric characters. If no string in the command, "H x" will read the stored header x.
Z	Same as pressing Zero Key.
T	Same as pressing Tare Key.
xT	Establish a preset Tare value in displayed unit. x = preset tare value. Sending 0T clears tare (if allowed).
PT	Prints Tare weight stored in memory.
PM	Print current application mode (weighing mode).
xM	Set current application mode to x. x depends on applications 1 – Weigh 2 – Count 3 – Percent 4 – Density
M	Scroll to the next enabled mode.
PU	Print Current weighing unit: g, kg, lb, oz, etc....
xU	Set balance to unit x: g, kg etc. 1 – g 2 – kg 3 – mg 4 – ct ...
U	Scroll to the next enabled unit.
ON	Brings out of Standby
OFF	Goes to Standby.
C	Begin Span Calibration
IC	Begin Internal Calibration.
AC	Abort Calibration.
PSN	Print Serial Number.
PV	Print Version: print name, software revision and LFT ON (if LFT is set ON).
x#	Set Counting APW (x) in grams. (must have APW stored)
P#	Print Counting application APW.
x%	Set Percent application reference weight (x) in grams. (must have reference weight stored)
P%	Print Percent application reference weight.
PTIME	Print current time.
PDATE	Print current date.
xTIME	Set Time x format: hh:mm:ss
xDATE	Set Date x format: mm/dd/yyyy
xS	0 = print unstable data, 1 = print stable only
xRL	Enable or disable OK response to non-print commands: x=0 to disable, x=1 to enable.
xT	Pre-tare the container weight (x) in grams.

## ACCESSORIES

For a complete listing of OHAUS printers and other accessories, contact OHAUS Corporation or visit [www.ohaus.com](http://www.ohaus.com).

## COMPLIANCE

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## DISPOSAL

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related. For disposal instructions in Europe, refer to [www.OHAUS.com/weee](http://www.OHAUS.com/weee). Thank you for your contribution to environmental protection.

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
Menu aggiunto	Funzione
Comunicazione  --RS232  --Tasso Baud	Impostare il tasso baud (bits per secondo).  <p>The screenshot shows a 'Baud Rate' menu with a grid of options: 1200, 2400, 4800, 9600 in the first row; 19200, 38400, 57600, 115200 in the second row. Below the grid are 'Back' and 'Exit' buttons.</p>
--Trasmissione	Impostare i data bits, stop bit, e parità.  <p>The screenshot shows a 'Transmission' menu with a grid of options: 7E1, 7E2, 7N1, 7N2 in the first row; 7O1, 7O2, 8N1, 8N2 in the second row. Below the grid are 'Back' and 'Exit' buttons.</p>
--Handshake	Impostare il metodo controllo flusso.  <p>The screenshot shows an 'RS-232' menu with three options: Baud Rate (with a crossed-out icon), Transmission (with a crossed-out icon), and Handshake (with a handshake icon). Below are 'None' and 'Xon/Xoff' buttons.</p>
-- Stampare le impostazioni	Per più informazioni, si prega di fare riferimento al Manuale d'Istruzioni STX.  <p>The screenshot shows a 'Communication' menu with two options: RS-232 (with a printer icon) and Print Settings (with a printer icon). Below are 'Back' and 'Exit' buttons.</p>

**Per Bilance PJX:**

Menu Navigazione	
<b>PRINT</b>	<b>RS232</b>
Header	Baud Rate
Date and Time	Transmission
Balance ID	Handshake
Balance Name	
User Name	
Project Name	
Application Name	
Result	
Gross Weight	
Net Weight	
Tare Weight	
Line Feed	

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<p><b>PRINT</b>                  Reset: <b>no, yes</b>                  Stable: <b>off, on</b>                  A.Print: <b>off</b>                            On.Stable                            interval                            accept                            continuous                  Content: Result (-&gt; <b>off, on</b>)                            Gross (-&gt; <b>off, on</b>)                            Net (-&gt; <b>off, on</b>)                            Tare (-&gt; <b>off, on</b>)                            Header (-&gt; <b>off, on</b>)                            Footer (-&gt; <b>off, on</b>)                            Mode (-&gt; <b>off, on</b>)                            Unit (-&gt; <b>off, on</b>)                            Info (-&gt; <b>off, on</b>)                  Layout: Format (-&gt;Single,<b>Multi</b>)                            Feed (-&gt;<b>Line</b>, 4 Lines, Form)                  End Print: Exit menu</p>	<p><b>RS232</b>                  Reset: <b>no, yes</b>                  Baud Rate: 1200...<b>9600</b>...115200                  Parity: 7 even, 7 odd, 7 none,                            <b>8 none</b>                  Stop bit: 1, 2                  Handshake: <b>none</b>, XOn-XOff                  Alternate command: Print (A...P...Z)                                            Tare (A...T...Z)                                            Zero (A...Z)                  End RS232: Exit menu</p>
---	---

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<b>메뉴 탐색</b>	
<b>PRINT</b> Header Date and Time Balance ID Balance Name User Name Project Name Application Name Result Gross Weight Net Weight Tare Weight Line Feed	<b>RS232</b> Baud Rate Transmission Handshake

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P/N 30269073 D © 2018 OHAUS Corporation, all rights reserved / todos los derechos reservados/ tous droits réservés/ Alle Rechte vorbehalten / Tutti i diritti riservati /版权所有 /모든 권리 보유 / Все права защищены /すべての権利予約