




pipetman®

Two Minute Inspection



1 STEP 1 Check the Records

- ▶ Use the serial number to identify the pipette and to determine its age
- ▶ Check laboratory records for the last service date


NOTE

The updated line of PIPETMAN Classic is identified with an underlined serial number starting with **QG**.



Letter	Year
A	1984 / 2006
B	1985 / 2007
C	1986 / 2008
D	1987 / 2009
E	1988 / 2010
G	1989 / 2011
H	1990 / 2012
J	1991 / 2013
K	1992 / 2014
L	1993 / 2015
M	1994 / 2016
N	1995 / 2017
P	1996 / 2018
Q	1997 / 2019
R	1998 / 2020
S	1999 / 2021
T	2000 / 2022
U	2001 / 2023
W	2002 / 2024
X	2003 / 2025
Y	2004 / 2026
Z	2005 / 2027

Letter	Month
A	January
B	February
C	March
D	April
E	May
G	June
H	July
J	August
K	September
L	October
M	November
N	December



2 STEP 2 General Appearance

CHECKED POINT	POSSIBLE CAUSES
Operating rod - bent? - corroded?	<ul style="list-style-type: none"> • Dropped • Lengthy immersion in corrosive liquid for decontamination • Lengthy exposure to corrosive vapors
Volumeter - dial alignment? - clarity of numbers?	<ul style="list-style-type: none"> • Autoclaving changed the appearance and function (the body must not be autoclaved)
Tip ejector - corroded? - broken?	<ul style="list-style-type: none"> • Lengthy immersion in corrosive liquid for decontamination • Lengthy exposure to corrosive vapors
Tip holder - physical or chemical damage?	<ul style="list-style-type: none"> • Repeated blows • Lengthy immersion in corrosive liquid for decontamination • Lengthy exposure to corrosive vapors

3 STEP 3 Check Functions

PROCEDURE	POSSIBLE CAUSES
Large volume adjustment 1 Set volume at maximum (i.e., nominal volume) assessing the movement of the friction ring 2 Activate the push button to test movement during aspirate and dispense strokes	<ul style="list-style-type: none"> • Irregular movement hitching, due to damage to the friction ring • No displacement bent operating rod • Jerky movement corroded, dirty, or scratched piston
Volumeter adjustment 1 Go through the entire range. The settings should correspond to the pipette's useful volume range (minimum to nominal volume)	<ul style="list-style-type: none"> • No adjustment autoclaving • Incorrect volume setting Misindexing; pipette adjustment screw has been incorrectly reassembled
Tip ejection system 1 Fit tip and depress tip ejector button 2 Observe function of tip ejector 3 Disassemble tip ejector	<ul style="list-style-type: none"> • No movement broken return spring • Improper fit not tight enough • Can't disassemble corrosion

4 STEP 4 Leak Test

PROCEDURE	POSSIBLE CAUSES
1 Fit PIPETMAN® DIAMOND tip 2 Set volume at maximum (i.e., nominal volume) 3 Pre-rinse by aspirating and dispensing water several times 4 Aspirate water 5 Hold the pipette in the vertical position for 20 seconds 6 For P2 to P200, re-immerses in the test liquid; fluid level in tip should remain constant 7 Observe if a drop or a leak appears at the orifice of the tip	<ul style="list-style-type: none"> • End of tip holder may be scratched/damaged (mechanical or chemical) • Improper fit • Use of non Gilson tips • Organic solvent, vapor pressure
Check these Tip holder - leak ? Tip - leak ?	

5 STEP 5 Disassembly - Reassembly

DISASSEMBLY	REASSEMBLY
1 Eject the tip 2 For PIPETMAN® Classic models manufactured before June 2019, with a blue tip ejector head, pull the tip ejector down For the updated models, with an underlined serial number starting from QG and a white tip ejector head, push laterally the tip ejector to slide and remove it 3 Unscrew the connecting nut 4 Separate the handle from the bottom part 5 Remove the piston from the tip holder	To avoid losing or damaging fragile parts, reassemble the pipette immediately Be sure to respect the correct order of parts: the piston seal should always be positioned before the O-ring
Check these Piston surface corroded, scratched, or, damaged Piston Seal and O-ring damaged (mechanical damage or chemical attack)	<p>NOTE You should never disassemble the body (handle) of the pipette.</p>