

Formaldehyde Sample Buffer

For best results, use within 1 year from date of manufacture.

Recommended procedure:

1. Thaw RNA samples on ice.
2. Thaw sample buffer at room temp.
3. Vortex sample buffer to mix thoroughly.
4. Dilute RNA sample 1:1 with Formaldehyde Sample Buffer.
5. Mix thoroughly.
6. Heat to 65°C for 5-15 minutes.
7. Cool for 1 minute on ice.
8. Load samples on gel.

Store: -20°C

For Research Use Only