

A 1 litre bottle of DISTEL DOSING SYSTEM concentrate solution generates 2000 measured doses of surface disinfectant. With the DISTEL DOSING SYSTEM every dose is the same volume as the previous one.

The DISTEL DOSING SYSTEM is intended for the cleaning and disinfection of non-invasive medical devices and hard surfaces within the hospital environment.

The DISTEL DOSING SYSTEM is effective against bacteria, yeast and viruses including coronaviruses from 3 minutes.

It is also destroys DNA and RNA.



ELIMINATE SINGLE-USE PLASTIC WIPES




The DISTEL DOSING SYSTEM is CE marked as a Class IIa Medical Device. *Dose onto a new wipe. For full instructions please see the DISTEL DOSING SYSTEM User Guide.

Designed to be repurposed and reused, the DISTEL DOSING SYSTEM bottles keep the environment in mind.

A 1 litre bottle of DISTEL DOSING SYSTEM concentrate delivers 10 x 1 litre bottles of working solution which equates to 200 x 5ml doses per bottle.

The DISTEL DOSING SYSTEM can be used with any wipe to create freshly made surface disinfectant wipes on the go. Choose a sustainably sourced wipe to reduce single-use plastic waste.





USE THE WIPE OF YOUR CHOICE

MANUFACTURED

CORONAVIRUSES

THE DISTEL DOSING SYSTEM IS INTENDED FOR THE CLEANING AND **DISINFECTION OF NON-INVASIVE MEDICAL DEVICES AND HARD** SURFACES WITHIN THE HOSPITAL ENVIRONMENT SUCH AS:



1x 500ml Bottle of Concentrate - 5x1Litre Dosing Bottles 6 x 5ml Dosing Heads - 1 x Working Solution Bottle Label

cache

BROUGHT TO YOU BY TRISTEL - T 01638 721500 - E mail@tristel.com - W www.tristel.com

stel patent information please visit: http://w our-patents.info/tristel - Copyright © Tristel Solutions - Mkt-Bro-1160-3 - April 2020