

Hybridisation Accessories

Hybridisation Incubators

Selectable speed rocking platform

Designed to cater for membrane-bound and slide-bound hybridisations it sits neatly in the base of the unit while still allowing a number of hybridisation tubes to be placed above. Made from stainless steel, the rocking platform enables membranes to be processed using a wave motion. Speed range of 0 to 60 oscillations per minute.

Membrane separators


Strong, re-usable porous mesh sheets for placing between membranes. Ensures even hybridisation when processing multiple hybridisation membranes in one tube. Available in packs of 5 (20 x 20cm)

Hybridisation tube rack

Useful storage facility for tubes during membrane loading or when not in use. Carries up to 3 large and 3 small tubes.

Tube holders

Tube holder with clips designed to carry 8 x 50ml tubes. Ideal for higher throughput of smaller samples (Hybrigene only).



Tube Holder



Rocking Platform

Ordering Information

Product Code	Description
FHB1/PLAT	Rocking platform 0-60opm (HB-1D only)
FHB4/PLAT	Rocking platform 0-60opm (Hybrigene only)
FMEM2020	Membrane separators 20 x 20cm (pack of 5)
FHTRACK	Tube rack holder (holds 3 large and 3 small tubes)
F50ML4TH	Tube carrier with clips to hold 8 x 50ml tubes