


Part No. 8050H-500

Description Disposable LiquiThief – 500mm

Material of Construction:

Tip	HDPE (virgin, FDA approved – CFR 177.1520)
Outer Tube	HDPE (virgin, FDA approved – CFR 177.1520)
Inner Rod	HDPE (virgin, FDA approved – CFR 177.1520)
Clutch	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Housing	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Seal	Silicone Rubber (virgin, FDA approved – CFR 177.2600)
Handle	HDPE (virgin, FDA approved – CFR 177.1520)



Method of Construction: Ultrasonic welding

Nominal Length

Fully closed:	523mm
Fully extended:	881mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: 1.5mm

Nominal Weight of Sampler: 83g

Assembly & Packing Environment: Class 100,000 Medical Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of Samplers per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:

Signature:



Name: Jonathan Fenner

Position: Director

Date: 29-Nov-12